

PLAC[®] Test Reagent Kit

Roche Cobas C 501 Chemistry Analyzer Application Sheet



This procedure is to be used in conjunction with the package insert for the PLAC Test Reagent Kit. Read and follow the package insert instructions and the Specimen Handling Best Practices carefully. Failure to follow the instructions may result in inaccurate results. Reagent barcode labels are not provided. For ordering information, contact Customer Service at 1-877-752-2837 or visit www.plactest.com. For further information, see the analyzer operator's manual.

PROCEDURE:

1. Program a user defined test with the parameters listed below. For more detailed instructions refer to the Roche Cobas C 501 User Instructions.
2. Transfer PLAC Test reagents into the appropriate reagent cartridge. Use Roche supplied evaporation caps to provide best on-board stability.
3. Load the cartridge as described in the User Operating Instructions.
4. Calibrate with PLAC calibrator set (**REF** 90108 or **REF** 10-0108).

NOTE:

- Do not report samples with results greater than 360 ng/mL.
- Do not dilute and retest.
- See PLAC Test turbidimetric immunoassay method package insert for further information. Refer to the sections: *Warnings & Precautions, Procedural Notes and Limitations.*

USER DEFINED PARAMETERS:

Analyze	Calib.		Range			Other	
Assay / Time / Point	2Point End		10	38	70	0	0
Wavelength (2nd / Pril.)	800		570				
Sample Volume			Cassette Configuration				
Norm.	4.8	0.0	0	Code			
Dec.	1.6	0.0	0	Expiration Days			
Inc.	4.8	0.0	0	*1			
Dilution			Reagent Volume				
o Water			R1	148	0	Inactive	
o Diluent			R2	0	0	Inactive	
			R3	48	0	Inactive	
Linearity Limit		%		%			
Prozone Limit						Inside	
Abs. Limit	32000		Increase				
Cell Detergent	*2		Stirring Level		8		
Stirring Setting			M1	M2	M3		
UP	Stirring		LOW	Stirring	Stirring	Stirring	

*1: User Defined

*2: Alkaline detergent

PLAC[®] Test Reagent Kit

Roche Cobas C 501 Chemistry Analyzer Application Sheet



Analyze	Calib.	Range	Other
Calibration Type <input type="text" value="Spline"/> Point <input type="text" value="5"/> Span <input type="text" value="5"/> Weight <input type="text" value="0"/> Update Type <input type="text" value="None"/> <input type="text" value="0"/> <input type="text" value="0"/>	SD Limit <input type="text" value="50"/> Duplicate Limit <input type="text" value="99"/> % <input type="text" value="32000"/> Abs. Sensitivity Limit <input type="text" value="-99999"/> <input type="text" value="99999"/> S1 Abs. Limit <input type="text" value="-32000"/> <input type="text" value="32000"/> <input type="checkbox"/> Auto Masking	Auto Calibration <input type="radio"/> Timeout Cassette <input type="text" value="Cancel"/> <input type="text" value="0"/> Day Changeover Cassette <input type="text" value="Cancel"/>	<input type="radio"/> QC Violation Method <input type="text" value="Blank"/> Rule <input type="text" value="1s"/> Control 1 <input type="text" value="None"/> Control 2 <input type="text" value="None"/> Control 3 <input type="text" value="None"/>

Analyze	Calib.	Range	Other																								
Application Code ### Unit ng/mL Report Name <input type="text" value="PLAC"/> Data Mode <input type="text" value="Active"/> <input type="checkbox"/> Automatic Return Technical Limit <input type="text" value="-9999"/> <input type="text" value="500"/> Report Limit <input type="text" value="-9999"/> <input type="text" value="99999"/> <input type="checkbox"/> Control Interval Time <input type="text" value="0"/> <input type="checkbox"/> Automatic QC On Board Stability <input type="text" value="1"/>	<input type="checkbox"/> Qualitative (1) <input type="text" value="0"/> (2) <input type="text" value="0"/> (3) <input type="text" value="0"/> (4) <input type="text" value="0"/> (5) <input type="text" value="0"/> (6) <input type="text" value="0"/>	Expected Values Male <table border="1"> <tr> <td></td> <td></td> <td><input type="text" value="-99999"/></td> <td><input type="text" value="999999"/></td> </tr> <tr> <td>99</td> <td>Year</td> <td><input type="text" value="-99999"/></td> <td><input type="text" value="999999"/></td> </tr> <tr> <td>100</td> <td>Year</td> <td><input type="text" value="-99999"/></td> <td><input type="text" value="999999"/></td> </tr> </table> Female <table border="1"> <tr> <td></td> <td></td> <td><input type="text" value="-99999"/></td> <td><input type="text" value="999999"/></td> </tr> <tr> <td>99</td> <td>Year</td> <td><input type="text" value="-99999"/></td> <td><input type="text" value="999999"/></td> </tr> <tr> <td>100</td> <td>Year</td> <td><input type="text" value="-99999"/></td> <td><input type="text" value="999999"/></td> </tr> </table>			<input type="text" value="-99999"/>	<input type="text" value="999999"/>	99	Year	<input type="text" value="-99999"/>	<input type="text" value="999999"/>	100	Year	<input type="text" value="-99999"/>	<input type="text" value="999999"/>			<input type="text" value="-99999"/>	<input type="text" value="999999"/>	99	Year	<input type="text" value="-99999"/>	<input type="text" value="999999"/>	100	Year	<input type="text" value="-99999"/>	<input type="text" value="999999"/>	Default Sex <input type="checkbox"/> Male <input type="checkbox"/> Female Range <input type="checkbox"/> Range 1 <input type="checkbox"/> Range 2 <input type="checkbox"/> Range 3
		<input type="text" value="-99999"/>	<input type="text" value="999999"/>																								
99	Year	<input type="text" value="-99999"/>	<input type="text" value="999999"/>																								
100	Year	<input type="text" value="-99999"/>	<input type="text" value="999999"/>																								
		<input type="text" value="-99999"/>	<input type="text" value="999999"/>																								
99	Year	<input type="text" value="-99999"/>	<input type="text" value="999999"/>																								
100	Year	<input type="text" value="-99999"/>	<input type="text" value="999999"/>																								

PLAC[®] Test Reagent Kit

Roche Cobas C 501 Chemistry Analyzer Application Sheet



Analyze	Calib.			Range		Other
Standards	(1)	(2)	(3)	(4)	(5)	(6)
Calibrator Code	*1	*1	*1	*1	*1	
Concentration	0.0	50.0	100.0	250.0	500.0	
Rack No.-Pos	*1	*1	*1	*1	*1	
Sample Volume	4.8	4.8	4.8	4.8	4.8	
Diluted S. Volume	0.0	0.0	0.0	0.0	0.0	0.0
Diluent Volume	0	0	0	0	0	0

*1: User Defined. Initial default parameters have been entered where appropriate.

PLAC[®] Test Reagent Kit

Roche Cobas C 501 Chemistry Analyzer Application Sheet



PERFORMANCE CHARACTERISTICS

Performance characteristics were established using the Roche Cobas C 501 chemistry analyzer.

On-board Analyzer Reagent Storage

Open bottles of reagents stored in the refrigerated compartment on the Cobas C501 analyzer should be stable up to 2 weeks. Use the reagent evaporation caps provided by Roche to give the best on-board stability. Laboratories should verify on-board stability on their own analyzers, under typical laboratory conditions.

Sensitivity

The clinical sensitivity of the assay is 12 ng/mL as determined by the limit of quantitation (the lowest concentration with acceptable precision).

The analytical sensitivity of the assay is 4 ng/mL, as calculated by interpolation of the mean plus two standard deviations of 20 replicates of the 0 ng/mL Lp-PLA₂ calibrator from the standard curve.

Assay Precision

Intra-assay and inter-assay variability were determined by testing two human serum samples with Lp-PLA₂ concentrations distributed throughout the calibration range of the assay. The data are summarized below:

Sample	Intra-assay		Inter-assay	
	Mean Concentration Lp-PLA ₂ (ng/mL)	% CV (n=20)	Mean Concentration Lp-PLA ₂ (ng/mL)	% CV (n=20)
Low Serum	189	3.1	150	6.0
High Serum	398	1.5	327	2.4

Linearity

From three pairs of serum samples with known high or low Lp-PLA₂ levels, a dilution series was prepared for each pair combination to assess linearity. Percent recoveries of the combined samples were determined as the measured value divided by the expected value, multiplied by 100. The average recovery was 100%, demonstrating linearity of the diluted samples over a range of 87 to 473 ng/mL Lp-PLA₂.

Interfering Substances

Endogenous substances found in blood were evaluated for interference in the assay. Four individual serum samples with Lp-PLA₂ values ranging from 276 to 441 ng/mL were spiked with potential interferents. No appreciable interference was observed for the following substances at the spiked levels tested.

Endogenous	
Potential Interferent	Test Concentration
Bilirubin	20 mg/dL
Cholesterol	500 mg/dL
Hemoglobin	500 mg/dL
Triglycerides	3000 mg/dL
Total Albumin*	~9000 mg/dL
* 5 g/dL albumin added to plasma pool of presumptively 4 g/dL albumin	

Method Comparison

Correlation studies were performed comparing the PLAC Test turbidimetric immunoassay to the values established on the Hitachi 917 and Olympus AU-400. Human serum samples with Lp-PLA₂ concentrations ranging from 117 to 481 ng/mL were tested. Results from linear regression analysis are shown below.

Linear Regression

Slope = 1.03

y-intercept = 18.5

Correlation coefficient r = 0.98

Number of samples = 30